

## New products

Filtration · Rapid Tests · Water Analysis · Chromatography · Bioanalysis  
Filtration · Schnellteste · Wasseranalytik · Chromatographie · Bioanalytik

### Filtration with suction filter Ashless filter MN 640 w – elevated border



<b>Grade:</b>	MN 640 w – elevated border
<b>Material:</b>	Refined pulp and linters, ashless (ash content <0,01%)
<b>Thickness:</b>	0,2 mm
<b>Filtration speed:</b>	9 s
<b>Basis weight:</b>	85 g/m <sup>2</sup>
<b>REF:</b>	202 901 (90 mm) / 202 902 (125 mm) / 202 903 (110 mm)
<b>Price:</b>	38,60 € / 46,60 € / 54,30 € per pack (valid until 31.12.2011)

**Content:** Pack with 2 x 50 filters

**Application:** Filtration with suction filters of 90 / 110 / 125 mm diameter

**Key customers:**

- Sewage plants
- Analytical laboratories
- Lab-suppliers

**Benefits:**

- The elevated border prevents residues to settle on the walls of the suction filter. Rinsing the suction filter is not necessary. The filter paper is taken from the suction filter and is available for further analysis right away.

**General information:** Filters with elevated border are a new innovation from MACHEREY-NAGEL and have been designed because of specific customer requests. They are ideal for filtration with suction filters. Using a classic filter circle in a suction filter leads to residues on the walls of the suction filter. Those residues have to be rinsed of the wall carefully so that they get onto the filter. Not rinsing the walls may lead to false results as not all residues are used for further analysis. MACHEREY-NAGEL filters with elevated border provide the solution to such problems. With the new design, all residues remain safely on the filter. Thus time-consuming and elaborate rinsing of suction filter walls is no necessary anymore. Therefore, the new filters effectively prevent false results.

Three different diameters are available:

- 90, 110 and 125 mm

If you are interested in the filters with a different paper grade, please contact us directly, as we could arrange a special production.

[www.mn-net.com](http://www.mn-net.com)



**MACHEREY-NAGEL GmbH & Co. KG** · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany  
**Germany and international:**  
 Tel.: +49 24 21 969-0  
 Fax: +49 24 21 969-199  
 E-mail: info@mn-net.com

**Switzerland:**  
**MACHEREY-NAGEL AG**  
 Tel.: +41 62 388 55 00  
 Fax: +41 62 388 55 05  
 E-mail: sales-ch@mn-net.com

**France:**  
**MACHEREY-NAGEL EURL**  
 Tel.: +33 388 68 22 68  
 Fax: +33 388 51 76 88  
 E-mail: sales-fr@mn-net.com

**USA:**  
**MACHEREY-NAGEL Inc.**  
 Tel.: +1 484 821 0984  
 Fax: +1 484 821 1272  
 E-mail: sales-us@mn-net.com

# New products

Filtration · Rapid Tests · Water Analysis · Chromatography · Bioanalysis  
Filtration · Schnellteste · Wasseranalytik · Chromatographie · Bioanalytik

**Accessories:** No accessories are necessary

If you require any further information please contact:

Technical support and customer service:

+49-(0)-2421 / 969-161  
+49-(0)-2421 / 969-174  
+49-(0)-2421 / 969-187  
e-mail: [csc@mn-net.com](mailto:csc@mn-net.com)

For orders:

+49 (0) 2421 / 969 0  
e-mail: [sales-de@mn-net.com](mailto:sales-de@mn-net.com)

